# Single Use Cryo Probe

INSTRUCTIONS FOR USE

|    |       |    |    |    |    |    |
|----|-------|----|----|----|----|----|
| EN | CN    | CS | DA | DE | ES | FR |
| IT | JP    | KO | LT | NL | NO | PL |
| PT | PT_BR | RO | SV |    |    |    |

 Keeler

|       |                                                                             |     |
|-------|-----------------------------------------------------------------------------|-----|
| EN    | Single Use Cryo Probe<br>Instructions For Use .....                         | 1   |
| CN    | 凯乐一次性冷冻探针<br>使用说明 .....                                                     | 8   |
| CS    | Kryosonda na jedno použití<br>Návod k použití .....                         | 15  |
| DA    | Kryoprobe til engangsbrug<br>Brugsanvisning.....                            | 22  |
| DE    | Kryosonde für den Einmalgebrauch<br>Gebrauchsanleitung.....                 | 29  |
| ES    | Criosonda monouso<br>Instrucciones De Uso.....                              | 36  |
| FR    | Sonde cryogénique à usage unique<br>Mode D'emploi .....                     | 43  |
| IT    | Sonda monouso Cryo<br>Istruzioni Per L'uso .....                            | 50  |
| JP    | キーラー単回使用クライオ・プローブ<br>取扱説明書.....                                             | 57  |
| KO    | 1회용 냉동탐침<br>사용설명서 .....                                                     | 64  |
| LT    | Vienkartinis kriozondas<br>Naudojimo instrukcija.....                       | 71  |
| NL    | Cryosonde voor eenmalig gebruik<br>Gebruiksaanwijzing .....                 | 78  |
| NO    | Kryoprobe til engangsbruk<br>Bruksanvisning.....                            | 85  |
| PL    | Jednorazowa sonda kriochirurgiczna<br>Instrukcja obsługi .....              | 92  |
| PT    | Sonda de Criocirurgia de Utilização Única<br>Instruções De Utilização ..... | 99  |
| PT_BR | Sonda Crio de Uso Único<br>Instruções de uso.....                           | 106 |
| RO    | Sondă criogenică de unică folosință<br>Instrucțiuni de utilizare .....      | 113 |
| SV    | Kryosond för engångsbruk<br>Bruksanvisning.....                             | 120 |

---

## CONTENTS

|                                                                                           |          |
|-------------------------------------------------------------------------------------------|----------|
| <b>1. INDICATIONS FOR USE .....</b>                                                       | <b>3</b> |
| <b>2. SAFETY .....</b>                                                                    | <b>3</b> |
| 2.1 WARNINGS AND CAUTIONS .....                                                           | 3        |
| 2.2 ATTENTION: IMPORTANT SAFETY INFORMATION .....                                         | 5        |
| 2.3 CONTRAINDICATION .....                                                                | 5        |
| <b>3. METHOD OF USE WITH KEELER DISPOSABLE PROBE ADAPTOR FOR MKI<br/>    CONSOLE.....</b> | <b>6</b> |
| <b>4. METHOD OF USE WITH KEELER CRYOMATIC MKII OR MKI CONSOLE .....</b>                   | <b>6</b> |
| <b>5. ENVIRONMENTAL CONDITIONS: .....</b>                                                 | <b>7</b> |
| <b>6. PACKAGING AND DISPOSAL INFORMATION.....</b>                                         | <b>7</b> |

|                           |                                                                              |                      |                                                                         |
|---------------------------|------------------------------------------------------------------------------|----------------------|-------------------------------------------------------------------------|
|                           | Consult instructions for use                                                 |                      | General warning sign                                                    |
|                           | Date of manufacture                                                          |                      | Warning: Electricity                                                    |
|                           | Manufacturer's name and address                                              |                      | Warning: Floor level obstacle                                           |
|                           | Country of manufacture                                                       |                      | Warning: Pressurized cylinder                                           |
|                           | Waste Electrical and Electronic Equipment (WEEE) recycling                   |                      | Warning; Low temperature/ Freezing conditions                           |
|                           | Keep dry                                                                     |                      | Batch Code                                                              |
|                           | Use-by date                                                                  |                      | Do not use if package is damaged                                        |
|                           | This way up                                                                  |                      | Fragile                                                                 |
|                           | Temperature limit                                                            |                      | Type BF applied part                                                    |
| <b>UK<br/>CA<br/>0120</b> | United Kingdom Conformity Assessed, with the Notified Body number for SGS UK | <b>CE<br/>1639</b>   | Conformité Européenne, with the Notified Body number for SGS Belgium NV |
| <b>EC</b> <b>REP</b>      | Authorised representative in the European Community                          | <b>CH</b> <b>REP</b> | Authorised representative in Switzerland                                |
|                           | Humidity limitation                                                          |                      | Atmospheric pressure limitation                                         |
| <b>REF</b>                | Catalogue number                                                             |                      | Do not re-use                                                           |
| <b>MD</b>                 | Medical device                                                               |                      | Sterilized using ethylene oxide                                         |
|                           | Translation                                                                  |                      |                                                                         |

**For use only with the Keeler Disposable Probe Adapter (2508-P-8026) and Cryomatic Console (2509-P-1000 or 2509-P-1000-JPN), or the Cryomatic MKII Console (2509-P-1010).**

**REF 2509-P-8033**

The Keeler Single Use Cryo Probe are designed and build in conformity with Directive 93/42/EEC, Regulation (EU) 2017/745 and ISO 13485 Medical Devices Quality Management Systems.

Classification: CE / UKCA: Class IIb

FDA: Class II

The information contained within this manual must not be reproduced in whole or part without the manufacturer's prior written approval. As part of our policy for continued product development we the manufacturer reserve the right to make changes to specifications and other information contained in this document without prior notice.

This IFU is also available on the Keeler UK and Keeler USA websites.

Copyright © Keeler Limited 2023. Published in the UK 2023.

## 1. INDICATIONS FOR USE

### Product Description

The Keeler Single Use Cryo Probe comprises a stainless steel stem with a plastic handle, connected to a long flexible plastic hose, comprising a high pressure tube within an outer protective tube, both of which are terminated within a plastic connector designed for insertion into the Keeler Cryomatic (with Disposable Probe Adaptor) and Keeler Cryomatic MKII.

**CAUTION: Federal Law restricts this device to sale by or on the order of a physician or practitioner.**

### Intended use

Keeler Single Use Cryo probes are for use in ophthalmic surgery in cryopexy for retinal detachment, cyclo destructive procedures in refractory glaucoma, cryo destruction of lash follicles for trichiasis and treatment of retinopathy of prematurity (ROP).

## 2. SAFETY

### 2.1 WARNINGS AND CAUTIONS

Please note that the proper and safe functioning of our instruments is only guaranteed if both the instruments and their accessories are exclusively from Keeler Ltd. The use of other accessories may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may lead to incorrect operation.

Observe the following precautions in order to ensure safe operation of the instruments.

### WARNINGS AND CAUTIONS

- Check your Single Use Cryo Probe for signs of transport / storage damage prior to use.
- Do not use in the presence of flammable gases / liquids, or in an oxygen rich environment.
- US Federal Law restricts this device to sale by or on the order of a physician or practitioner.
- Only use a Keeler Disposable Probe Adapter (2508-P-8026) or the Cryomatic MKII Console (2509-P-1010).
- This device is intended to be used only by suitably trained and authorised healthcare professionals experienced in the use of ophthalmic surgical equipment.
- Both Single Use Cryo Probe and Disposable Probe Adapter instructions are to be strictly followed to avoid any usage errors.
- This is a sterile device for single use on one patient only. For one use only. Re-use of a single device may result in cross-contamination and infection.
- Do not put this device in any sterilising / decontaminating solution.
- Do not use after the validity date indicated after the symbol on the label. Use of the device after the expiration date could lead to the risk of infection.
- Fragile. Avoid knocks and handle with care.
- This product should not be immersed in fluid.
- No modification to this equipment is allowed.
- Using the device out of specification or when there is a fault can lead to electric shock, overheating or burns (hot or cold).
- Do not use this device if the packaging is in any way broken or damaged. If packaging/sterile barrier is damaged, risk of infection can occur.
- Do not repair and use a damaged product.
- Usage of this device for Ophthalmic surgery(s) using cryotherapy (dependant on the type of surgery being performed) may lead to the following complications/side effects such as transient conjunctival edema, eyelid edema, meibomian gland dysfunction, skin depigmentation, recurrence of the deformation of the eyelashes, scar progression and rarely, retinal tear and detachment, visual loss, lid notching, acceleration of symblepharon formation, xerosis, cellulitis, conjunctival hyperaemia, conjunctival chemosis, macular drag, macula pucker, retinal folds, conjunctival lacerations, vitreous haemorrhages/retinal haemorrhages, increase in IOP caused by gas injection during pneumatic retinopexy, damage to surrounding tissue area, dispersion of RPE cells after cryopexy, diplopia, and cryo-induced chorioretinal.
- In the event of failure to defrost the user may be unable to remove the probe tip from the eye.
- Only the areas of the eye to be treated should make contact with the probe tip or probe shaft.
- Do not use if the hose is twisted or kinked.

- Test the Single Use Cryo Probe before use.
- Do not use probes which are leaking.
- Do not obstruct the vent holes on the Single Use probe (located on the tube).
- Do not attempt to eject the probe while the footswitch is depressed.
- After use treat as contaminated waste and dispose of according to local hospital regulations.
- Do not re-sterilise.
- When the tip is not in view, use the orientation pip on the handle to identify the tip position.
- Use only genuine Keeler approved parts and accessories or device safety and performance may be compromised.
- Keep out of the reach of children.
- For indoor use only. This device should be used in a sterile hospital operating theatre environment indicated after the symbol on the label.

## 2.2 ATTENTION: IMPORTANT SAFETY INFORMATION

1. Test the probe operation prior to use.
2. Connect the Cryo Probe and depress the foot pedal for 20-30 seconds. The probe should freeze within 5 seconds of pressing the pedal and defrost within 5 seconds of releasing it.

**If the above test fails, or you notice a leak, or other malfunction, turn off the gas supply and follow precisely, the instruction on how to disconnect the Single Use Cryo Probe/ console, as described in the relevant instruction for use. Return the faulty Single Use Cryo Probe to Keeler for analysis and exchange.**

3. Take a new Single Use Cryo Probe and refer to the user instructions detailed above.

**If the probe tip fails to defrost within 5 seconds of the footswitch being released, immediately switch off the console, and apply sterile water to the tip contact area. This will defrost the tip within a few seconds. Return the faulty Single Use Cryo Probe to Keeler for analysis and exchange.**

4. Take a new Single Use Cryo Probe and refer to the user instructions detailed above.
5. Please ensure that the probe is decontaminated prior to return and is accompanied by a decontamination certificate indicating that it is safe to handle. It is vital that the LOT number of the probe is recorded on the accompanying documentation – this will enable us to identify potential problems with specific batches and take appropriate action.

## 2.3 CONTRAINDICATION

There is no restriction to patient population this device can be used with other than those outlined in the contraindications stated in Section 2.1 Warnings and Cautions.

### 3. METHOD OF USE WITH KEELER DISPOSABLE PROBE ADAPTOR FOR MKI CONSOLE

1. Insert the Adapter for the Single Use Cryo Probe into the standard probe socket of the Cryo Unit. Leave or hold the hinged safety clip facing upwards.
2. Insert the Single Use Cryo Probe connector into the Adapter. Push it firmly all the way in until you hear the click of the sprung safety lock.






3. Lower the security clip to ensure that it secures the tubing at the end of the connector.
4. If running software version vG or earlier the Keeler Cryomatic will automatically purge for 90 seconds.
5. If running software version vH or later the Keeler Cryomatic will display the 'Run' symbol. 
6. Commence operation of Cryo Unit according to manufacturer's instructions.
7. Test Single Use Cryo Probe before applying to the eye.
8. At the end of the operation, carefully follow shut down procedure of the Cryo Unit. Vent the gas by purging according to the equipment instructions.
9. To remove the Single Use Cryo Probe, first press the sprung safety lock, then lift the hinged safety clip and remove the used probe assembly. Please ensure it is disposed of safely and is not re-used.
10. Use the enclosed peel-off labels to identify the product used on patient/facility records.

### 4. METHOD OF USE WITH KEELER CRYOMATIC MKII OR MKI CONSOLE

Refer to the Keeler Cryomatic MKII or MKI Cryosurgical System instructions.

## 5. ENVIRONMENTAL CONDITIONS:

| USE                  |                                                                                                                                                                                                                                                       |                       |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|
|                      |    |                       |
|                      | Shock (without packing)                                                                                                                                                                                                                               | 10 g, duration 6 ms   |
| STORAGE CONDITIONS   |                                                                                                                                                                                                                                                       |                       |
|                      |    |                       |
| TRANSPORT CONDITIONS |                                                                                                                                                                                                                                                       |                       |
|                      |    |                       |
|                      | Vibration, sinusoidal                                                                                                                                                                                                                                 | 10 Hz to 500 Hz: 0.5g |
|                      | Shock                                                                                                                                                                                                                                                 | 30 g, duration 6 ms   |
|                      | Bump                                                                                                                                                                                                                                                  | 10 g, duration 6 ms   |

**Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of your Member State.**

## 6. PACKAGING AND DISPOSAL INFORMATION

### Disposal of old electrical and electronic equipment

This symbol on the product or on its packaging and instructions indicates that this product shall not be treated as household waste.

To reduce the environmental impact of WEEE (Waste Electrical Electronic Equipment) and minimise the volume of WEEE entering landfills we encourage at product end of life that this equipment is recycled and reused.

**If you need more information on the collection reuse and recycling then please contact B2B Compliance on 01691 676124 (+44 1691 676124). (UK only).**

**Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of your Member State.**

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Copyright © Keeler Limited 2023. Publikováno v UK v roce 2023.
Copyright © Keeler Limited 2023. Publiceret 2023 i Storbritannien.
Copyright © Keeler Limited 2023. Veröffentlicht in GB 2023.
Copyright © Keeler Limited 2023. Publicado en el Reino Unido 2023.
Copyright © Keeler Limited 2023. Publié au Royaume Uni en 2023.
Copyright © Keeler Limited 2023. Pubblicato nel Regno Unito 2023.